#### Welcome to the




# Clinical diagnostics with automated SNP and CNV analysis using varvis™

Dr. Yvonne Schmitz, Product Management yvonne.schmitz@limbus-medtec.com


### What to expect

- Brief introduction to the varvis genomics platform
  - Components and workflow
  - Data regulation, protection and availability
- Examples: real case studies from routine diagnostics
  - Presented by Dr. Konrad Platzer, Human Genetics, University Hospital Leipzig
- Hands-on practical session
  - Puzzle over some real case and play with the software yourself
- Solutions and discussion
- Final remarks






### varfeed: Variant calling pipeline for NGS




### CNV analysis from targeted panels




# Validated for clinical application

# Performed on 430 TruSight Cancer Panels on MiSeq Comparison to MLPA




Parameter	Symbol	Value
True positive	TP	117
True negative	TN	48236
False positive	FP	10
False negative	FN	0
Sensitivity/Recall	TPR=TP/(TP+FN)	100.00%
Specificity	SPC=TN/(TN+FP)	99.98%


#### CNV visualization








# Varvis Workflow








#### Public data: Annotation sources + Version


- 1000 Genomes: 2014\_10 Afrian, American, EastAsian, European, SouthAsian, Total Frequencies
- Clinvar: 201711: Accession Number, Significance, Status
- **ESP: 1.0** EVS Total Frequency
- ExAC: 03 African, American, EastAsian, European NonFinnish, Finnish, Other, SouthAsian, Total Frequencies
- GnomAD: 2.0.1 Female, HomInd, Male, Total GnomExs
- SCSNV: 1.1 ScSnvAda, ScSnvRf
- dbNSFP: 3.4 Insilico Prediction Rankings, including GERP++, Mutation Taster, Sift and many more
- dbSNP: 142

#### Public data: Annotation sources + Version

U25-A1 - Index



## allexes Protects Patient Data by Design



## Medical device as cloud platform

#### Data privacy

- BDSG
- EU Regulation 2016/679




Certified data center Data security assessment

#### Medical approval

- Regulation 93/42/EEC
- MPG
- ISO 13485
- ISO 14971

Quality management system Clinical validation Risk management





#### UNIVERSITÄT LEIPZIG

Medizinische Fakultät



#### variant evaluation in NGS diagnostics (clinical examples in live demo + Hands-on practical session without slides)

**Konrad Platzer, MD** 

Institute of Human Genetics University of Leipzig, Germany



# Final remarks

### Summary – key benefits of varvis™

#### C Up to date

Use continuously updated public data from sources like **gnomAD** and **ClinVar**.

#### **₩** 6000+

Access more than 6000 diligently curated WES cases including all variants and phenotype descriptions through our network.

#### Overnight express

No matter how many samples or how many sequencing runs you have: our fully automated process delivers results overnight.

#### 80% savings

Our clinically validated CNV analysis provides significant cost savings compared to ligation based CNV detection methods.

#### **C€** Approved

varvis™ is a medical device Class I according to MDD 93/42/EEC.

#### Deep phenotyping

varvis™ enables standardized phenotyping using ontologies like HPO to improve the accuracy of your diagnosis.

### Thank you for attending our workshop!

Do you have additional questions or want a live demo or varvis?

Are you interested in the clinical validation of your CNV data?

Would you like to test varvis with your own data?

Please visit <a href="www.varvis.com">www.varvis.com</a> for more information Or contact yvonne.schmitz@limbus-medtec.com